


# The SAFE Study

Safety and Efficacy Analysis of FRED Embolic Device in Aneurysm Treatment: One-year Clinical and Anatomical Results

*A single arm, multicenter, prospective, GCP study of 103 patients with intracranial aneurysms treated with the Flow-Redirection Intraluminal Device (FRED®) system at 13 French neurovascular centers.*

## Study Purpose

To evaluate the safety and efficacy of the Flow-Redirection Intraluminal device, (FRED®) and FRED® Jr. devices in real-world patients treated for intracranial aneurysms.



## Key Findings at 12 month follow-up

**73.3%**

Complete Occlusion

**81.1%**

Adequate Occlusion

**1.9%**

Mortality


\* No device related mortality

**2.9%**

Morbidity

## Progressive Occlusion Observed Over Time

\*Measured by 3-Point Raymond-Roy Scale



## Author's Conclusion

"SAFE study analysis at one year confirms the excellent safety profile of FRED® device for aneurysm treatment with low morbidity and mortality rates (2.9% and 1.9% respectively) and demonstrates its efficacy (adequate occlusion in 81.1%)."

## Safety Findings

**2.2%**

Retreatment

**95.1%**

Good Clinical  
Outcomes (mRs 0-2)

### SAFE STUDY POPULATION

Population	103 aneurysms / 103 patients
Age	25 to 80 years old (mean = 52.4)
<b>ANEURYSM LOCATION</b>	
Supraclinoid ICA	68.9% (71)
Cavernous ICA	14.6% (15)
MCA	7.8% (8)
ACA or ACom	8.7% (9)
<b>ANEURYSM SIZE</b>	
Small (<10mm)	68.9% (71)
Large (10-24mm)	28.2% (29)
Giant (>/=24mm)	2.9% (3)
Wide Neck	96.1% (99)
<b>TREATMENT CHARACTERISTICS</b>	
Technical Success	95.1%

### Safety Discussion:

"Noticeably, the safety of aneurysm treatment is now, at least with the FRED® device, very close to the safety of standard coiling as reported for example in ATENA (with one-month morbidity and mortality being 1.7% and 1.4% respectively)"

Study limitations are outlined in the full article, available at: <https://jn.is.bmj.com/content/early/2018/10/08/neurintsurg-2018-014261.share>

Pierot L, Spelle L, Berge J, et al. SAFE study (Safety and Efficacy Analysis of FRED Embolic Device in Aneurysm Treatment): One-year Clinical and Anatomical Results. J Neurointervent Surg Epub ahead of print. October 9, 2018. doi = 10.1136/neurint

RX Only: For health care professional use only. This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopy equipment. The FRED® flow diversion device should only be used by physicians who have received appropriate training for the device. The FRED® System is not approved in the U.S. by the FDA.

#### INDICATIONS FOR USE:

The system is intended for endovascular embolization of intracranial neurovascular aneurysms. The FRED® and FRED® Jr. system may also be used with embolic coils for the treatment of intracranial neurovascular lesions.

**CE 0297**



[www.microvention.com](http://www.microvention.com)

**MicroVention USA**  
35 Enterprise  
Aliso Viejo, CA 92656  
USA  
+1 714 247 8000  
PH 1 800 990 8368

**MicroVention United Kingdom**  
Suite 3, The Barracks Building  
10 Cliords Fort, North Shields  
Tyne and Wear, NE30 1JE  
United Kingdom  
PH +44 (0) 191 258 6777  
F +44 (0) 191 258 5999

**MicroVention Europe SARL**  
30 bis, rue du Vieil Abreuvoir  
78100 Saint-Germain-en-Laye  
France  
RCS Versailles B 440 775674  
00029  
PH +33 (1) 39 21 52 17  
F +33 (1) 39 21 16 01

**MicroVention  
Deutschland GmbH**  
Hildebrandtstr. 4 F  
D-40215 Düsseldorf  
Germany  
PH +49 211 210 798-0  
F +49 211 210 798-29