

PLATINUM COILS

MicroPlex® Coil System

COMPLEX COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Standard	18	Cosmos® 18	24	68
Standard	18	Cosmos® 18	22	63
Standard	18	Cosmos® 18	20	65
Standard	18	Cosmos® 18	18	59
Standard	18	Cosmos® 18	16	52
Standard	18	Cosmos® 18	14	51
Standard	18	Cosmos® 18	13	47
Standard	18	Cosmos® 18	12	43
Standard	18	Cosmos® 18	11	39
Standard	18	Cosmos® 18	10	36
Standard	18	Cosmos® 18	9	31
Standard	18	Cosmos® 18	8	27
Standard	18	Cosmos® 18	7	23
Standard	18	Cosmos® 18	6	19
Soft	10	Cosmos® 10	10	36
Soft	10	Cosmos® 10	9	33
Soft	10	Cosmos® 10	8	37
Soft	10	Cosmos® 10	8	25
Soft	10	Cosmos® 10	7	31
Soft	10	Cosmos® 10	7	22
Soft	10	Cosmos® 10	6	26
Soft	10	Cosmos® 10	6	18
Soft	10	Cosmos® 10	5	22
Soft	10	Cosmos® 10	5	15
Soft	10	Cosmos® 10	4	12
Soft	10	Cosmos® 10	4	8
Soft	10	Cosmos® 10	3	6
Soft	10	Cosmos® 10	2.5	4
Soft	10	Cosmos® 10	2	2
Extra Soft	10	HyperSoft® 3D	6	20
Extra Soft	10	HyperSoft® 3D	5	20
Extra Soft	10	HyperSoft® 3D	5	15
Extra Soft	10	HyperSoft® 3D	5	10
Extra Soft	10	HyperSoft® 3D	4	15
Extra Soft	10	HyperSoft® 3D	4	12
Extra Soft	10	HyperSoft® 3D	4	8
Extra Soft	10	HyperSoft® 3D	4	6
Extra Soft	10	HyperSoft® 3D	3.5	8
Extra Soft	10	HyperSoft® 3D	3.5	5
Extra Soft	10	HyperSoft® 3D	3	10
Extra Soft	10	HyperSoft® 3D	3	8
Extra Soft	10	HyperSoft® 3D	3	6
Extra Soft	10	HyperSoft® 3D	3	4
Extra Soft	10	HyperSoft® 3D	2.5	8
Extra Soft	10	HyperSoft® 3D	2.5	6
Extra Soft	10	HyperSoft® 3D	2.5	4
Extra Soft	10	HyperSoft® 3D	2	8
Extra Soft	10	HyperSoft® 3D	2	6
Extra Soft	10	HyperSoft® 3D	2	4
Extra Soft	10	HyperSoft® 3D	2	3
Extra Soft	10	HyperSoft® 3D	2	2
Extra Soft	10	HyperSoft® 3D	1.5	4
Extra Soft	10	HyperSoft® 3D	1.5	3
Extra Soft	10	HyperSoft® 3D	1.5	2
Extra Soft	10	HyperSoft® 3D	1	3
Extra Soft	10	HyperSoft® 3D	1	2

HELICAL COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Extra Soft	10	HyperSoft® Helical	6	8
Extra Soft	10	HyperSoft® Helical	5	8
Extra Soft	10	HyperSoft® Helical	4	8
Extra Soft	10	HyperSoft® Helical	4	6
Extra Soft	10	HyperSoft® Helical	3	8
Extra Soft	10	HyperSoft® Helical	3	6
Extra Soft	10	HyperSoft® Helical	3	4
Extra Soft	10	HyperSoft® Helical	2.5	6
Extra Soft	10	HyperSoft® Helical	2.5	4
Extra Soft	10	HyperSoft® Helical	2	6
Extra Soft	10	HyperSoft® Helical	2	4
Extra Soft	10	HyperSoft® Helical	2	3
Extra Soft	10	HyperSoft® Helical	2	2
Extra Soft	10	HyperSoft® Helical	1.5	4
Extra Soft	10	HyperSoft® Helical	1.5	3
Extra Soft	10	HyperSoft® Helical	1.5	2
Extra Soft	10	HyperSoft® Helical	1.5	1
Extra Soft	10	HyperSoft® Helical	1.5	4
Extra Soft	10	HyperSoft® Helical	1.5	3
Extra Soft	10	HyperSoft® Helical	1.5	2
Extra Soft	10	HyperSoft® Helical	1.5	1

SPECIALTY FILLING COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Soft	N/A	VFC®	15-20	60
Soft	N/A	VFC®	15-20	40
Soft	N/A	VFC®	10-15	40
Soft	N/A	VFC®	10-15	30
Soft	N/A	VFC®	6-10	30
Soft	N/A	VFC®	6-10	20
Extra Soft	N/A	VFC®	3-6	15
Extra Soft	N/A	VFC®	3-6	10
Extra Soft	N/A	VFC®	3-6	6
Extra Soft	N/A	VFC®	1-3	3

HYDROGEL COILS

HydroCoil® Embolic System

COMPLEX COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Standard	18	HydroFrame® 18	20	48
Standard	18	HydroFrame® 18	18	50
Standard	18	HydroFrame® 18	16	44
Standard	18	HydroFrame® 18	14	45
Standard	18	HydroFrame® 18	13	47
Standard	18	HydroFrame® 18	12	43
Standard	18	HydroFrame® 18	11	39
Standard	18	HydroFrame® 18	10	36
Standard	18	HydroFrame® 18	9	31
Standard	18	HydroFrame® 18	8	27
Standard	18	HydroFrame® 18	7	23
Standard	18	HydroFrame® 18	6	19
Soft	10	HydroFrame® 10	10	36
Soft	10	HydroFrame® 10	9	31
Soft	10	HydroFrame® 10	8	33
Soft	10	HydroFrame® 10	8	17
Soft	10	HydroFrame® 10	7	28
Soft	10	HydroFrame® 10	7	15
Soft	10	HydroFrame® 10	6	19
Soft	10	HydroFrame® 10	6	12
Soft	10	HydroFrame® 10	5	20
Soft	10	HydroFrame® 10	5	15
Soft	10	HydroFrame® 10	5	10
Soft	10	HydroFrame® 10	4	15
Soft	10	HydroFrame® 10	4	10
Soft	10	HydroFrame® 10	4	8
Soft	10	HydroFrame® 10	4	5
Extra Soft	10	HydroSoft® 3D	3	10
Extra Soft	10	HydroSoft® 3D	3	8
Extra Soft	10	HydroSoft® 3D	3	6
Extra Soft	10	HydroSoft® 3D	3	4
Extra Soft	10	HydroSoft® 3D	2.5	8
Extra Soft	10	HydroSoft® 3D	2.5	6
Extra Soft	10	HydroSoft® 3D	2.5	4
Extra Soft	10	HydroSoft® 3D	2	8
Extra Soft	10	HydroSoft® 3D	2	6
Extra Soft	10	HydroSoft® 3D	2	4
Extra Soft	10	HydroSoft® 3D	2	3
Extra Soft	10	HydroSoft® 3D	2	2
Extra Soft	10	HydroSoft® 3D	1.5	4
Extra Soft	10	HydroSoft® 3D	1.5	3
Extra Soft	10	HydroSoft® 3D	1.5	2
Extra Soft	10	HydroSoft® 3D	1	3
Extra Soft	10	HydroSoft® 3D	1	2

HELICAL COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Soft	10	HydroSoft® 10	6	20
Soft	10	HydroSoft® 10	6	15
Soft	10	HydroSoft® 10	6	10
Soft	10	HydroSoft® 10	6	8
Soft	10	HydroSoft® 10	5	20
Soft	10	HydroSoft® 10	5	15
Soft	10	HydroSoft® 10	5	10
Soft	10	HydroSoft® 10	5	8
Soft	10	HydroSoft® 10	5	6
Soft	10	HydroSoft® 10	4	10
Soft	10	HydroSoft® 10	4	8
Soft	10	HydroSoft® 10	4	6
Soft	10	HydroSoft® 10	4	4
Soft	10	HydroSoft® 10	3	10
Soft	10	HydroSoft® 10	3	8
Soft	10	HydroSoft® 10	3	6
Soft	10	HydroSoft® 10	3	4
Soft	10	HydroSoft® 10	2.5	6
Soft	10	HydroSoft® 10	2.5	4
Soft	10	HydroSoft® 10	2	8
Soft	10	HydroSoft® 10	2	6
Soft	10	HydroSoft® 10	2	4
Soft	10	HydroSoft® 10	2	3
Soft	10	HydroSoft® 10	2	2
Soft	10	HydroSoft® 10	2	1
Soft	10	HydroSoft® 10	1.5	4
Soft	10	HydroSoft® 10	1.5	2

SPECIALTY FILLING COILS				
Softness	Coil Size	Product Item Description	Diameter (mm)	Length (cm)
Soft	10	HydroFill®	12	30
Soft	10	HydroFill®	10	30
Soft	10	HydroFill®	9	30
Soft	10	HydroFill®	8	20
Soft	10	HydroFill®	7	20
Soft	10	HydroFill®	6	20
Soft	10	HydroFill®	6	15
Soft	10	HydroFill®	5	15
Soft	10	HydroFill®	5	10
Soft	10	HydroFill®	4	10
Soft	10	HydroFill®	4	8
Soft	10	HydroFill®	3	6
Soft	10	HydroFill®	2	6

Available on the V-Trak® Advanced Delivery System

INDICATIONS FOR USE:

The Hydrocoil® Embolic System (HES) and MicroPlex Coil System (MCS) are intended for the endovascular embolization of intracranial aneurysms and other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae. The HES and MCS are also intended for vascular occlusion of blood vessels within the neurovascular system to permanently obstruct blood flow to an aneurysm or other vascular malformation and for arterial and venous embolizations in the peripheral vasculature.

RX Only: Federal (USA) law restricts this device to sale by or on the order of a physician. This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopy equipment. The MicroPlex® Coil System and HydroCoil® Embolic System Coils should only be used by physicians who have received appropriate training for the device.

MICROVENTION, MicroPlex, Hydrocoil, HydroFrame, HydroFill, HydroSoft, Cosmos, HyperSoft, VFC, and V-Trak are registered trademarks of MicroVention, Inc. Refer to Instructions for Use, contraindications and warnings for additional information. Federal (USA) law restricts this device for sale by or on the order of a physician.

©2018 MicroVention, Inc. MM621(i) APLA 6/18